



### FEATURES

- Powered neck loop to connect to any Smartphone, iPhone, MP3, Kindle...that has a 3.5mm headset socket
- Rechargeable batteries with up to 2 weeks standby time
- Switch between music and phone calls at the touch of a button.
- Incoming voice amplification
- Gain up to 30 dB
- Built-in microphone
- In use indication by LED (Green)
- Low battery indicator by LED (Red)
- Charge indicator by LED (Blue)
- Simple connection: 3.5 mm
- Switch power ON/OFF
- Breakout connector

### TECHNICAL SPECIFICATION

- Micro USB charging connector
- USB cable 1m (included)
- Length of cable of iLoop+: 1.15m
- Operating time: 5 hours
- Powered by 2 AAA 1.2V/1000mA niMH (included)
- Size: 3.5 x 2.2 x 7 cm
- Compatible with:  
CL8360 and CL8500 and all other mobile phone with a 3.5mm jack